# **BRILLIANT** EverGlow<sup>®</sup> The glow of the art





## THE GLOW OF THE ART

#### THE DESIRE TO HAVE AN EASY-TO-PROCESS, AESTHETIC AND DURABLE FILLING MATERIAL IS SIGNIFICANT.

BRILLIANT EverGlow is a material that meets the demand without compromise. The state-of-the art composite is suitable for anterior and posterior restorations and offers the right mix of simple and practical handling, excellent blend-in properties and long-lasting brilliance.

- Exceptional polishability and long-lasting brilliance
- Aesthetic single-shade restorations
- Smooth consistency provides ideal handling
- Good wettability on the conditioned tooth surface



## EXCEPTIONAL POLISHABILITY AND LONG-LASTING BRILLIANCE

### BRILLIANT EVERGLOW IS A UNIVERSAL SUBMICRON HYBRID COMPOSITE.

Pre-polishing (i.e. with DIATECH COMPREPOL PLUS)

Due to its sophisticated filler composition, BRILLIANT EverGlow immediately shows an exceptionally smooth surface which has a satin shine directly after placing the filling. This makes polishing a minor issue and highly aesthetic restorations can be performed in next to no time. In addition, the versatile filling material excels through its pronounced gloss retention.




Gloss polish (i.e. with DIATECH COMPOSHINE PLUS)



Easy to polish to a glossy final result instantly

#### **GLOSS RETENTION**<sup>1</sup>

After exposure in the toothbrush simulator, measured in gloss units (GU)



23.1 % very good good modest bad

A field test with 39 dentists from Europe and the US has confirmed the outstanding polishability of BRILLIANT EverGlow: 100 % of the participants gave a positive judgment, 76.9 % of them giving the highest possible score.

Even after loading for 6,000 cycles in a toothbrush simulator, BRILLIANT EverGlow retains its initial gloss very well.


<sup>1</sup>Source: internal data | \* not a trademark of COLTENE

### **EVALUATION OF POLISHABILITY**<sup>1</sup>

## AESTHETIC SINGLE-SHADE RESTORATIONS WITH UNIVERSAL<sup>1</sup> SHADES

BRILLIANT EverGlow harmoniously integrates into existing surroundings, with the application of a single shade already providing unerring highly aesthetic restorations. An additional enamel layer can often be dispensed with. Exept in case of very young patients or when remodelling especially translucent incisal edges, can one resort to enamel shades. This was made possible by further developing the adaptive «Duo Shade» system of COLTENE. Due to its excellent optical blend-in properties, each shade covers two VITA\* shades at the same time, resulting in A1 / B1 or A2 / B2, for example.

\* VITA is not a trademark of Coltène/Whaledent



### CLINICAL CASES WITH UNIVERSAL<sup>1</sup> SHADES ONLY



## A FLEXIBLE SYSTEM FOR ALL INDICATIONS

BRILLIANT EverGlow is a sensible and truly versatile system giving the dentist freedom to achieve highly aesthetic restorations. In addition to the universal shades there are two translucentshades which provide the option of an individual characterisation of particularly translucent enamel areas. The BRILLIANT EverGlow shade system is finally completed with three opaque shades. With these strongly discoloured teeth can be masked, providing the basis for an aesthetic restoration.



- Simple shade management: A single shade provides two VITA shades ("Duo Shade" System)
- 7 universal composite shades, suitable for anterior and posterior regions
- Additional enamel and opaque shades for aesthetically more complex restorations



All pictures: Courtesy of the Styleitaliano Team <sup>1</sup>"Dentin" shades for single-shade restorations

### SMOOTH CONSISTENCY PROVIDES IDEAL HANDLING

Its smooth consistency and dimensional stability make BRILLIANT EverGlow a state-of-the-art filling material easy and convenient to apply in cavities. Elaborately modelled cusps, ledges and contact points remain intact and do not slump. High resistance of the composite to prolonged exposure to surgery lighting completes ideal handling. The smooth consistency and excellent wettability of BRILLIANT EverGlow make it ideal for luting COMPONEER, the Direct Composite Veneering System by COLTENE.




## GOOD WETTABILITY ON THE CONDITIONED TOOTH SURFACE

#### **STICKINESS TO THE INSTRUMENT<sup>1</sup>**


Measured in Newton

80

Due to its good wettability, BRILLIANT EverGlow offers excellent adhesion to bonded tooth surfaces rather than sticking to the instruments. This has also been confirmed in an international user survey: over 90 % of participants rated the wetting properties as being «very good» or «good».



WETTABILITY ON BOND<sup>1</sup>




## THE KEY TO BRILLIANT RESTORATIONS

## THE BRILLIANT EVERGLOW FILLER TECHNOLOGY

BRILLIANT EverGlow's filler technology is all geared to provide an ideal mixture of handling convenience, mechanical strength and long-lasting aesthetics. It is based on three core facts: submicron barium glass fillers care for excellent polishability and gloss retention where as BRILLIANT EverGlow's pre-polymerised fillers lower volume shrinkage and improve sculpting properties of the paste without impairing polishability and gloss retention thanks to their sophisticated composition ideally harmonizing with the composite. Above all, the filler surfaces are treated to ensure optimal compound strength within the resin matrix and support gloss retention. The result is a state-of-the-art composite that allows for "brilliant" restorations.


# SIMPLY IMPRESSIVE

SEM picture of BRILLIANT EverGlow



#### **MINIMAL ABRASION 1**

Relative loss of volume (dimensionless)



<sup>1</sup> Source: internal data | \* not a trademark of COLTENE

#### **MINIMAL ABRASION**

Compared to other conventional composites, BRILLIANT EverGlow displays significantly lower abrasion values (see graph above).

#### HIGH COMPRESSIVE STRENGTH

Ideal for restorations of posterior teeth due to high resistence to mastication forces.

#### LOW WATER ABSORPTION

The optimised composition leads to lower water absorption, which in turn leads to excellent long-term plour stability.

### ORDER INFORMATION



BRILLIANT EverGlow Kits	Opaque Shades			Universal Shades							Translucent Shades	
	Opaque Bleach (OBL)	Opaque A1 (OA1)	Opaque A3 (OA3)	Bleach (BL)	A1 / B1	A2/B2	A3 / D3	A3.5/B3	C2/C3	A4 / C4	Bleach Translucent (BL Trans)	Translucent (Trans)
Trial Kits												
6001 9720   4 Syringes à 3g						1	1	1				1
6001 9721   32 Tips à 0,2g						8	8	8				8
Anterior Kits												
6001 9724   4 Syringes à 3g				1	1						1	1
6001 9725   32 Tips à 0,2g				8	8						8	8
System Kits												
6001 9722   9 Syringes à 3g				1	1	1	1	1	1	1	1	1
6001 9723   60 Tips à 0,2g				5	8	8	8	8	5	5	5	8
Regular Kits												
6001 9726   9 Syringes à 3g		1	1		1	1	2	1	1			1
6001 9727   60 Tips à 0,2g		5	5		5	10	12	10	5			8
Anterior Kits Large												
6001 9728   9 Syringes à 3g	1	1		2	2	1					1	1
6001 9729 60 Tips à 0,2g	5	5		12	12	10					8	8

#### BRILLIANT EverGlow Refills Tips (20 × 0,2 g)

6001 9746	Opaque Bleach (OBL)		
6001 9747	Opaque A1 (OA1)		
6001 9748	Opaque A3 (OA3)		
6001 9706	Bleach (BL)		
6001 9707	A1/B1		
6001 9708	A2/B2		
6001 9709	A3/D3		
6001 9710	A3.5/B3		
6001 9716	C2/C3		
6001 9717	A4/C4		
6001 9718	Translucent (Trans)		
6001 9719	Bleach Translucent (BL Trans)		

#### BRILLIANT EverGlow Refills Syringes (1 × 3 g)

nemis synnges (r x s g/			
6001 9743	Opaque Bleach (OBL)		
6001 9744	Opaque A1 (OA1)		
6001 9745	Opaque A3 (OA3)		
6001 9697	Bleach (BL)		
6001 9698	A1/B1		
6001 9699	A2/B2		
6001 9700	A3/D3		
6001 9701	A3.5/B3		
6001 9702	C2/C3		
6001 9703	A4/C4		
6001 9704	Translucent (Trans)		
6001 9705	Bleach Translucent (BL Trans)		





© 2016 Coltène/Whaledent AG - www.coltene.com

Coltène/Whaledent Ltd. The President Suite Kendal House, Victoria Way Burgess Hill, West Sussex RH15 9NF / UK T +44 1444 235486 F +44 1444 870640 info.uk@coltene.com Coltène/Whaledent AG Feldwiesenstr. 20 9450 Altstätten / Switzerland T +41 71 757 53 00 F +41 71 757 53 01 info.ch@coltene.com Coltène/Whaledent Inc. 235 Ascot Parkway Cuyahoga Falls, Ohio 44223 / USA T +1 330 916 8800 F +1 330 916 7077 info.us@coltene.com

